

# PRODUCT INFORMATION Thermo Scientific Optimized DyNAzyme Buffer

#\_\_\_\_

Lot: \_

## Ordering information

Component	#F-511S	#F-511L
10X Optimized DyNAzyme Buffer	2 x 1.5 ml	4 x 1.5 ml

Store at -20°C

**|∥** 37

### www.thermoscientific.com/pcr

## Description

10X Thermo Scientific Optimized DyNAzyme Buffer should be diluted 1:10 in the final PCR mix. 1X buffer contains 10 mM Tris-HCI (pH 8.8 at 25°C),1.5 mM MgCl<sub>2</sub>, 50 mM KCI, 0.1% Triton X-100.

#### Technical support

US: techservice.genomics@thermofisher.com Europe, Asia, Rest of World: techservice.emea.genomics@thermofisher.com Web: www.thermoscientific.com/pcr

#### PRODUCT USE LIMITATION

This product has been developed and is sold exclusively for research purposes and in vitro use only. This product has not been tested for use in diagnostics or drug development, nor are they suitable for administration to humans or animals. © 2012 Thermo Fisher Scientific, Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

