

NativePAGE™ Sample Buffer (4X)

(Novex®)

货号 BN2003

包装 10 ml

NativePAGE™ Sample Buffer (4X) is a solution used to solubilize native protein samples prior to use with the NativePAGE™ Gel System. The buffer is used along with 10% DDM or 5% Digitonin detergents to solubilize hydrophobic or membrane proteins. This buffer is also included in the NativePAGE™ Sample Prep Kit. The samples prepared with 10% DDM, 5% Digitonin, or the NativePAGE™ Sample Prep Kit are compatible with NativePAGE™ Novex Bis-Tris Gels for native electrophoresis showing increased resolution and reduced streaking.

详细说明

常用规格

Gel Type:	Native Gel
Buffer Type:	Sample Loading Buffers
Concentrated:	4 X
Product Size:	10 ml
Gel Compatibility:	NativePAGE™ Gels
Shipping Condition:	Wet Ice
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

内容和储存

Store at 2°C to 8°C