

Recombinant Human FGF-9 (Fibroblast Growth Factor-9)

Publication Number MAN0004378 Rev. 1.00

| Catalog Number: | PHG0194 |
|-----------------------------------|---|
| Quantity: | 20 μg |
| Lot Number: | See product label. |
| Molecular Weight: | 23.4 kDa, 207 amino acid residues. |
| Purity: | ≥95% by SDS PAGE and HPLC analyses. |
| Biological Activity: | $ED_{50} \le 0.5 \text{ ng/mL}$ ($\ge 2 \times 10^6 \text{ units/mg}$), determined by the dose dependent stimulation of thymidine uptake by BaF3 cells expressing FGF-receptors. FGF–9 is a heparin binding growth factor that stimulates the proliferation of glial cells, as well as other cells expressing the FGF–9 receptor. Cellular sources of FGF-9 include prostatic stromal cells, neurons, astrocytes, oligodendrocytes, fetal skeletal myoblasts, ventricular myocardium, gut lumenal epithelium and olfactory bulb. FGF–9 is known alternatively as glial-activating factor (GAF). |
| Formulation: | Lyophilized from 1xTris Buffer, pH 7.8. Filtered through a 0.2 micron filter prior to lyophilization. |
| Endotoxin: | <0.1 ng/μg of FGF–9 |
| Production: | Produced in E. coli and purified by sequential chromatography. |
| Reconstitution Recommendation: | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. This lyophilized protein should be reconstituted in sterile, distilled water to a concentration of 0.1–1 mg/mL. <i>Do not vortex</i> . This solution can be stored at 2°C to 8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing a carrier protein, such as 0.1% BSA, and store in working aliquots at –20°C to -80°C. |
| Suggested Working Dilutions: | The optimal concentration should be determined for each specific application. |
| Storage: | The lyophilized protein is stable at room temperature for 1 month, but should be kept at –20°C for long term storage, preferably desiccated. Working aliquots stored with a carrier protein are stable for at least 3 months at –20°C to –80°C. Avoid repeated freeze/thaw cycles. |
| Expiration Date: | See product label. |
| References: | Giri, D. et al. (1999) J. Cell Physiol. 180:53. Nakamura, S. et al. (1999) Glia 28:53. Todo, T. et al. (1998) Brain Res. 783:179. Colvin, J.S. et al. (1999) Dev. Dyn. 216:72. |

Explanation of Symbols

The symbols present on the product label are explained below:

| Symbol | Description |
|------------|---|
| REF | Catalog Number |
| RUO | Research Use Only |
| | Use by |
| *** | Manufacturer |
| [-] | Without, does not contain |
| from Light | Protect from light |
| <u> </u> | Directs the user to consult instructions for use (IFU), accompanying the product. |

| Symbol | Description |
|----------|--|
| LOT | Batch code |
| IVD | In vitro diagnostic medical device |
| X | Temperature limitation |
| EC REP | European Community authorized representative |
| [+] | With, contains |
| <u> </u> | Consult accompanying documents |

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

 $\label{lem:continuous} \textbf{For Research Use Only. Caution: Not for human or animal the rapeutic or diagnostic use.}$

 $Manufacturing \ site: 7335 \ Executive \ Way \ | \ Frederick, \ MD \ 21704 \ | \ Toll \ Free \ in \ USA \ 800.955.6288$

© 2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit www.lifetechnologies.com/support or email techsupport@lifetech.com www.lifetechnologies.com

Revision Date 20 May 2011

