

Recombinant Human Bone Morphogenetic Protein-7, Active (BMP-7, Active)

Publication Number MAN0003498 Rev. 2.00

Catalog Number:	PHC7204	
Quantity:	10 µg	
Lot Number:	See product label.	
Molecular Weight:	28.8 kDa homodimer, each subunit contains 116 amino acid residues (corresponding to amino acid residues 316 431 of the full-length BMP–7 precursor)	
Purity:	≥98% by SDS-PAGE and HPLC analyses.	
Biological Activity:	$ED_{50} = 0.02-0.04~\mu g/mL$, determined by the dose dependent induction of alkaline phosphatase production by murine ATDC5 cells.	
Formulation:	Lyophilized from 10 mM sodium citrate. Sterile filtered through a 0.2 micron filter prior to lyophilization.	
Endotoxin:	<0.1 ng/µg	
Production:	Produced in CHO cells and purified by sequential chromatography.	
Reconstitution Recommendation:	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in sterile, distilled water, to a concentration of 0.1–1.0 mg/mL. For extended storage, we recomment to further dilute in a buffer containing carrier protein (ex: 0.1% BSA) and store in working aliquots at –20°C to -80°C.	
Suggested Working Dilutions:	The optimal concentration should be determined for each specific application.	
Storage:	This lyophilized preparation is stable at room temperature for up to 1 month, but should be kept at –20°C for long term storage, preferably desiccated. Upon reconstitution, the preparation is stable for up to one week at 2°C to 8°C Working aliquots stored with a carrier protein are stable for at lease 7 months at –20°C to –80°C. Avoid repeated freeze/thaw cycles.	
Expiration Date:	See product label.	
References:	 Dudley, A.T. and E.J. Robertson (1997) Overlapping expression domains of bone morphogenetic protein family members potentially account for limited tissue defects in BMP7 deficient embryos. Dev. Dyn. 208:349–362. Dudley, A. T., et al. (1995) A requirement for bone morphogenetic protein-7 during development of the mammalian kidney and eye. Genes Dev. 9:2795–2807. Luo, G., et al. (1995) BMP-7 is an inducer of nephrogenesis, and is also required for eye development and skeletal patterning. Genes Dev. 9:2808–2820. Hu, M.C., et al. (2004) p38MAPK acts in the BMP7-dependent stimulatory pathway during epithelial cell morphogenesis and is regulated by Smad1. J. Biol. Chem. 279(13):12051–12059. 	

Explanation of Symbols

The symbols present on the product label are explained below:

Symbol	Description
REF	Catalog Number
RUO	Research Use Only
	Use by
***	Manufacturer
[-]	Without, does not contain
from Light	Protect from light
<u> </u>	Directs the user to consult instructions for use (IFU), accompanying the product.

Symbol	Description
LOT	Batch code
IVD	In vitro diagnostic medical device
1	Temperature limitation
EC REP	European Community authorized representative
[+]	With, contains
<u> </u>	Consult accompanying documents

Limited Use Label License: Research Use Only

The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product only to perform internal research for the sole benefit of the purchaser. No right to resell this product or any of its components is conveyed expressly, by implication, or by estoppel. This product is for internal research purposes only and is not for use in commercial applications of any kind, including, without limitation, quality control and commercial services such as reporting the results of purchaser's activities for a fee or other form of consideration. For information on obtaining additional rights, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

For Research Use Only. Caution: Not for human or animal therapeutic or diagnostic use.

Manufacturing site: 7335 Executive Way | Frederick, MD 21704 | Toll Free in USA 800.955.6288

© 2011 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

For support visit www.lifetechnologies.com/support or email techsupport@lifetech.com www.lifetechnologies.com



