Certificate of Analysis

pGL4.74[hRluc/TK] Vector:

 Part No.
 Size

 E692A
 20μg



Instructions for use of this product can be found in the pGL4 Vectors Technical Manual #TM259, available online at: **www.promega.com/protocols**

Description: The pGL4.74[hRluc/TK] Vector(a-d) encodes the luciferase reporter gene hRluc (Renilla reniformis) and is designed for high expression and reduced anomalous transcription. The pGL4 Vectors are engineered with fewer consensus regulatory sequences and a synthetic gene, which has been codon optimized for mammalian expression.

The pGL4.74[hRluc/TK] Vector contains the hRluc luciferase reporter gene and an HSV-TK promoter and can be used as an expression control or a co-reporter vector.

Concentration: 1µg/µl.

GenBank® Accession Number: AY738230.

Storage Buffer: The pGL4.74[hRluc/TK] Vector is supplied in 10mM Tris-HCI (pH 7.4), 1mM EDTA.

Storage Conditions: See the product information label for storage temperature recommendations and expiration date information

Usage Notes: Concentration gradients may form in frozen products and should be dispersed upon thawing. Mix well prior to use.

Quality Control Assays

Contaminant Assays

Contaminating Nucleic Acids: RNA, single-stranded DNA and chromosomal DNA are not evident in specified quantities of the vector as determined by agarose gel electrophoresis.

Nuclease Assay: Following incubation of 1µg of the vector in Restriction Enzyme Buffer at 37°C for 16–24 hours, no evidence of nuclease activity is detected by agarose gel electrophoresis.

Physical Purity: $A_{260}/A_{280} \ge 1.80$, $A_{260}/A_{250} \ge 1.05$.

Functional Assays

Identity Assay: The vector has been sequenced completely and has 100% identity with the published sequence available at: www.promega.com/vectors/

Restriction Digestion: The functional purity of the vector DNA is verified by successful digestion with restriction enzymes at the optimal temperature for one hour. Samples are examined by agarose gel electrophoresis, comparing cut and uncut vector DNA with marker DNA.

Signed by:

J. Stevens, Quality Assurance

(a) BY USE OF THIS PRODUCT, RESEARCHER AGREES TO BE BOUND BY THE TERMS OF THIS LIMITED USE LABEL LICENSE. If the researcher is not willing to accept the terms of this label license, and the product is unused, Promega will accept return of the unused product and provide the researcher with a full refund.

Stevens

Researchers may use this product for research use only, no commercial use is allowed. "Commercial use" means any and all uses of this product and derivatives by a party for money or other consideration and may include but is not limited to use in: (1) product manufacture; and (2) to provide a service, information or data; and/or resale of the product or its derivatives, whether or not such product or derivatives are resold for use in research. Researchers shall have no right to modify or otherwise create variations of the nucleotide sequence of the luciferase gene except that researchers may: (1) create fused gene sequences provided that the coding sequence of the resulting luciferase gene except that researchers may: (1) create fused gene sequences provided that the coding sequence of the resulting luciferase gene has no more than four deoxynucleotides missing at the affected terminus compared to the intact luciferase gene sequence, and (2) insert and remove nucleic acid sequences in splicing research predicated on the inactivation or reconstitution of the luminescence of the encoded luciferase. No other use or transfer of this product or derivatives is authorized without the prior express written consent of Promega. In addition, researchers must either: (1) use luminescence activity of this product and its derivatives; or (2) contact Promega to obtain a license for use of the product and its derivatives. Researchers may transfer derivatives to there for research use provided that at the time of transfer a copy of this label license is given to the recipients and recipients agree to be bound by the terms of this label license. With respect to any uses outside this label license, including any diagnostic, therapeutic or prophylactic uses, please contact Promega for supply and licensing information. PROMEGA MAKES NO REPRESENTATIONS OR WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING FOR MER-CHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WITH REGARDS TO THE PRODUCT. The terms of this label license

(b)Patent pending.

(c)U.S. Pat. No. 7,906,282 and European Pat. No. 1341808

Part# 9PIE692 Revised 8/13





Promega Corporation	
2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Toll Free	800-356-9526
Fax	608-277-2516
Internet	www.promega.com

PRODUCT USE LIMITATIONS, WARRANTY, DISCLAIMER

Promega manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Promega products contain chemicals which may be armful if misused. Due care should be exercised with all Promega products to prevent direct human contact.

Each Promega product is shipped with documentation stating specifications and other technical information. Promega products are warranted to meet or exceed the stated specifications. Promega's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Promega makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO PROMEGA PRODUCTS. In no event shall Promega be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Promega products to perform in accordance with the stated specifications.

© 2004–2013 Promega Corporation. All Rights Reserved.

GenBank is a registered trademark of the US Department of Helath and Human Services. All specifications are subject to change without prior

All specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

Part# 9PIE692 Printed in USA. Revised 8/13



Features List and Map for the pGL4.74 [hRluc/TK] Vector

The following features are present in the vector based on nucleotide sequence.

HSV-TK promoter	27-779
hRluc reporter gene	815–1750
1 0	
SV40 late poly(A) signal	1784–2005
Reporter Vector primer 4 binding region	2071–2090
CoIEI-derived plasmid replication origin	2330
Synthetic β-lactamase (Ampr) coding region	3119-3979
Synthetic poly(A) signal/transcriptional pause site	4084-4237
Reporter Vector primer 3 binding region	4186-4205




Figure 1. pGL4.74 [hRluc/TK] Vector circle map and sequence reference points.