# Technical Data Sheet

# **AEC Substrate Set**

### **Product Information**

**Material Number:** 551951 Size: 10 Plate(s) 51-2577KC Component: **Description:** AEC Substrate 100 mL (1 ea) Size:

Component: 51-2578KC **Description:** AEC Chromogen 3 mL (1 ea) Size:

## Description

3-Amino-9-ethylcarbazole (AEC) is a substrate for peroxidase and is a widely-used chromogen for ELISPOT.

# Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

## **Application Notes**

#### Application

ELISPOT Routinely Tested

#### Recommended Assay Procedure:

No more than 15 minutes prior to use, mix one drop (20 µL) of AEC Chromogen with each 1 mL of AEC Substrate, vortex or mix. Discard any remaining prepared AEC Substrate after use.

Warning: AEC Substrate (MN 51-2577KC) contains 7.99% methanol (w/w).

## Hazard statements

Combustible liquid.

Harmful if swallowed.

May be harmful in contact with skin.

May cause damage to organs.

## Precautionary statements

Wear protective gloves / eye protection.

Wear protective clothing.

Avoid breathing mist/vapours/spray.

IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.

IF ON SKIN: Wash with plenty of water.

Danger: AEC Chromogen (MN 51-2578KC) contains 97.8% N,N-dimethylformamide (w/w).

# Hazard statements

Flammable liquid and vapour.

Harmful if inhaled.

Causes serious eye irritation.

May damage fertility or the unborn child.

# Precautionary statements

Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

Wear protective clothing / eye protection.

Wear protective gloves.

IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Dispose of contents/container in accordance with local/regional/national/international regulations.

## **BD Biosciences**

bdbiosciences.com

551951 Rev. 3

**United States Asia Pacific** 866.979.9408 32.2.400.98.95 0120.8555.90 65.6861.0633 877.232.8995 55.11.5185.9995

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD



## **Suggested Companion Products**

Catalog Number	Name	Size	Clone
557630	HRP Streptavidin for ELISPOT	1 mL	(none)

## **Product Notices**

1. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.



bdbiosciences.com

 
 United States
 Canada
 Europe
 Japan

 877.232.8995
 866.979.9408
 32.2.400.98.95
 0120.8555.90
 **Asia Pacific** Latin America/Co 65.6861.0633 55.11.5185.9995 Latin America/Caribbean

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD

