## **Technical Data Sheet**

# PE Mouse Anti-Human IL-4

## **Product Information**

Material Number:	554516
Size:	0.1 mg
Concentration:	0.2 mg/ml
Clone:	8D4-8
Immunogen:	Recombinant Human IL-4
Isotype:	Mouse IgG1, ĸ
Reactivity:	QC Testing: Human
Storage Buffer:	Aqueous buffered solution containing $\leq 0.09\%$ sodium azide.

## Description

The 8D4-8 monoclonal antibody reacts with human interleukin-4 (IL-4). The immunogen used to raise the 8D4-8 hybridoma was recombinant human IL-4. The 8D4-8 antibody binds to an epitope that is different than the epitope recognized by the MP4-25D2 antibody (Cat. No. 554485).

Clone 8D4-8 displays an increased amount of non-specific binding to dead cells when compared to the clone MP4-25D2. It is recommended to use a fixable viability dye in conjunction with this clone.



Expression of IL-4 by stimulated human peripheral blood mononuclear cells (PBMC). Human PBMC were stimulated with immobilized anti-human CD3 antibody (10  $\mu$ g/ml for coating culture wells; Cat. No. 555329), soluble anti-human CD28 (2  $\mu$ g/ml final concentration; Cat. No. 555725) recombinant human IL-2 (10 ng/ml final concentration; Cat. No. 554603) and recombinant human IL-4 (20 ng/ml final concentration; Cat. No. 555726) recombinant human IL-2 (210 ng/ml final concentration; Cat. No. 554605) for 2 days. The cells were subsequently cultured in medium containing recombinant human IL-2 and recombinant human IL-4 for 3 days. Finally, the cells were harvested and stimulated for 4 hours with PMA (Sigma, Cat. #P-8139) and ionomycin (Sigma, Cat. #C-9275) in the presence of GolgiStop<sup>TM</sup> (2  $\mu$ M final concentration; Cat. No. 554724). The cells were harvested, stained with FITC-anti- CD4 (Cat. No. 555346), fixed, permeabilized, and subsequently stained with 0.015  $\mu$ g of PE-mouse anti-human IL-4 antibody (PE-8D4-8, Cat. No. 554516) by using Pharmingen's staining protocol (see image, left panel). To demonstrate specificity of staining, the binding of PE-8D4-8 antibody was blocked by preincubation of the antibody with recombinant human IL-4 (Cat. No. 554605; middle panel) or by preincubation of the fixed/permeabilized cells with unlabelled 8D4-8 antibody (10  $\mu$ g, Cat. No. 554515; right panel) prior to staining. The quadrant markers for the bivariate dot plot were set based on the autofluorescence controls and verified using the recombinant cytokine and unlabelled antibody blocking controls.

## Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with R-PE under optimum conditions, and unconjugated antibody and free PE were removed.

## **Application Notes**

#### Application

Intracellular staining (flow cytometry)	Routinely Tested	
---	------------------	--

**Recommended Assay Procedure:** 

**Immunofluorescent Staining and Flow Cytometric Analysis:** The PE-conjugated 8D4-8 antibody can be used for multicolor immunofluorescent staining and flow cytometric analysis to identify and enumerate IL-4-producing cells within mixed cell populations (see image, left panel). For optimal immunofluorescent staining for flow cytometric analysis, titration of the antibody is strongly recommended. For specific methodology, please visit our web site, www.bdbiosciences.com, and go to the protocols section or the chapter on intracellular staining in the Immune Function Handbook. A suitable mouse IgG1 isotype control for assessing the level of background staining on paraformaldehyde - fixed/saponin -permeabilized human cells is PE-MOPC-21 (Cat. No. 554680); use at comparable concentrations to antibody of interest.

#### **BD Biosciences**

bdbiosciences.com							
United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean		
877.232.8995	800.268.5430	32.2.400.98.95	0120.8555.90	65.6861.0633	55.11.5185.9995		
For country cor	For country contact information, visit <b>bdbiosciences.com/contact</b>						
Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD							



## Suggested Companion Products

Catalog Number Name		Size	Clone	
554680	PE Mouse IgG1, κ Isotype Control	0.1 mg	MOPC-21	
555329	Purified NA/LE Mouse Anti-Human CD3	0.5 mg	UCHT1	
555725	Purified NA/LE Mouse Anti-Human CD28	0.5 mg	CD28.2	
554715	BD Cytofix/Cytoperm Plus Kit (with BD GolgiStop)	250 tests	(none)	
554605	Recombinant Human IL-4	5 µg	(none)	
554603	Recombinant Human IL-2	10 µg	(none)	
555062	HiCK-2 Human Cytokine Positive Control Cells	1.0 ml	(none)	

## **Product Notices**

Since applications vary, each investigator should titrate the reagent to obtain optimal results. 1.

- 2. An isotype control should be used at the same concentration as the antibody of interest.
- 3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at 4. www.bdbiosciences.com/colors.
- 5. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

#### References

Bird C, Wadhwa M, Thorpe R. Development of immunoassays for human interleukin 3 and interleukin 4, some of which discriminate between different recombinant DNA-derived molecules. Cytokine. 1991; 3(6):562-567. (Clone-specific)

Prussin C, Metcalfe DD. Detection of intracytoplasmic cytokine using flow cytometry and directly conjugated anti-cytokine antibodies. J Immunol Methods. 1995; 188(1):117-128. (Methodology)

#### **BD Biosciences**

bdbiosciences.com United States Canada

Asia Pacific 65.6861.0633 
 Canada
 Europe
 Japan

 800.268.5430
 32.2.400.98.95
 0120.8555.90
877.232.8995 For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2014 BD

Latin America/Caribbean

55.11.5185.9995

