# **Technical Data Sheet**

# Alexa Fluor® 488 Mouse anti-PLK1 (pT210)

#### **Product Information**

 Material Number:
 558446

 Size:
 50 tests

 Vol. per Test:
 20 μl

 Clone:
 K50-483


Immunogen: Phosphorylated Human PLK1 Peptide

Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

# Storage Buffer: **Description**

Polo-like kinase (PLK1) is a serine/threonine kinase with structural similarities to Drosophila's Polo kinase and the Cdc5p of Saccharomyces cerevisiae. Like its invertebrate counterparts, PLK1 activity is required for DNA synthesis and is regulated throughout the cell cycle. Furthermore, PLK1 is highly expressed in primary tumors. It associates with the mitotic spindle during mitosis suggesting that, in addition to its role during S phase, PLK1 may play a role during chromosome segregation. This is consistent with its potential role in cancer development. Threonine 210 (T210) is one of the major phoshorylation sites in activated PLK1 obtained from human mitotic cells.

The K50-483 monoclonal antibody recognizes the phosphorylated T210 of human PLK1.



Analysis of PLK1 (pT210) in human epithelioid carcinoma. Hela S3 cells (ATCC CCL 2.2) were either stimulated with Nocodazole at 37°C for 12-16 hours (shaded histogram) or unstimulated (open histogram). The cells were fixed (BD Cytofix™ buffer, Cat. No. 554655) for 10 minutes at 37°C, then permeabilized (BD Phosflow™ Perm Buffer III, Cat. No. 558050) on ice for at least 30 minutes, and then stained with Alexa Fluor® 488 Mouse anti-PLK1 (pT210). Simultaneous staining with 7-AAD (Cat. No. 559925) demonstrates that the cells with upregulated PLK1 phosphorylation are in the G2/M phase of the cell cycle (data not shown). Flow cytometry was performed on a BD FACSCalibur™ flow cytometry system.

#### **Preparation and Storage**

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated to Alexa Fluor® 488 under optimum conditions, and unreacted Alexa Fluor® 488 was removed.

## **Application Notes**

Application

Intracellular staining (flow cytometry)  Routinely Tested	Intracellular staining (flow cytometry)	Routinely Tested	
---	---	------------------	--

## **Suggested Companion Products**

Catalog Number	Name	Size	Clone
558050	Perm Buffer III	125 ml	(none)
554655	Fixation Buffer	100 ml	(none)

#### **BD Biosciences**

bdbiosciences.com

 United States
 Canada
 Europe
 Japan
 Asia Pacific
 Latin America/Caribbean

 877.232.8995
 800.979.9408
 32.53.720.550
 0120.8555.90
 65.6861.0633
 55.11.5185.9995

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be constructed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be help responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only, Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD



558446 Rev. 2 Page 1 of 2

#### **Product Notices**

- This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1 × 10<sup>6</sup> cells in a 100-μl experimental sample (a test).
- 2. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 3. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 4. Alexa Fluor® is a registered trademark of Molecular Probes, Inc., Eugene, OR.
- 5. Alexa Fluor® 488 fluorochrome emission is collected at the same instrument settings as for fluorescein isothiocyanate (FITC).
- 6. The Alexa Fluor®, Pacific Blue™, and Cascade Blue® dye antibody conjugates in this product are sold under license from Molecular Probes, Inc. for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430), Pacific Blue™ dye, and Cascade Blue® dye are covered by pending and issued patents.
- 7. For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at www.bdbiosciences.com/colors.
- 8. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

#### References

Ando K, Ozaki T, Yamamoto H, et al. Polo-like kinase 1 (Plk1) inhibits p53 function by physical interaction and phosphorylation. *J Biol Chem.* 2004; 279(24):25549-25561. (Biology)

Jang Y-J, Ma S, Terada Y, Erikson RL. Phosphorylation of threonine 210 and the role of serine 137 in the regulation of mammalian Polo-like kinase. *J Biol Chem.* 2002; 277(46):44115-44120. (Biology)

Peter B, Gleixner K, Cerny-Reiterer S, et al. Polo-like kinase-1 as a novel target in neoplastic mast cells: demonstration of growth-inhibitory effects of small interfering RNA and the polo-like kinase-1 targeting drug Bl 2536. *Haematologica*. 2011; 96(5):672-680. (Clone-specific: Immunocytochemistry (cytospins)) Yamashiro S, Yamakita Y, Totsukawa G, et al. Myosin phosphatase targeting subunit1 regulates mitosis by antagonizing polo-like kinase1. *Dev Cell*. 2008; 14(5):787-797. (Clone-specific: Immunofluorescence, Western blot)

#### **BD Biosciences**

bdbiosciences.com

 United States
 Canada
 Europe
 Japan
 Asia Pacific
 Latin America/Caribbean

 877.232.8995
 800.979.9408
 32.53.720.550
 0120.8555.90
 65.6861.0633
 55.11.5185.9995

For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be constructed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be help responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD



558446 Rev. 2