

43-4423 2 ml HRP-Avidin

Lot No. See product label Exp. Date: See product label

PEROXIDASE-AVIDIN (FROM EGG WHITE)

FORM:

2 ml liquid. The product is in 10 mM phosphate buffer-0.15 M NaCl, pH 7.4, containing enzyme conjugate stabilizer, and 0.1% proclin as preservative. The protein concentration of peroxidase conjugated avidin is 2.5 mg/ml.

AVIDIN (from egg white) is conjugated with highly purified horseradish peroxidase (RZ > 3.0)

RECOMMENDED WORKING DILUTION FOR ELISA:

1:2,000 to 1:5,000. Working dilutions for specific applications should be determined by the investigator to obtain the best conditions.

STORAGE:

Store at 2-8°C.

Explanation of symbols				
Symbol	Description	Symbol	Description	
REF	Catalogue Number	LOT	Batch code	
RUO	Research Use Only	IVD	In vitro diagnostic medical device	
X	Use by	ł	Temperature limitation	
***	Manufacturer	EC REP	European Community authorised representative	
[-]	Without, does not contain	[+]	With, contains	
Street Light	Protect from light	Æ	Consult accompanying documents	
$\prod i$	Directs the user to consult instructions for use (IFU), accompanying the product.			

For research use only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

PI434423 (Rev 03/10) DCC-10-0838

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.



43-4423 2 ml HRP-Avidin

Lot No. See product label Exp. Date: See product label

PEROXIDASE-AVIDIN (FROM EGG WHITE)

FORM:

2 ml liquid. The product is in 10 mM phosphate buffer-0.15 M NaCl, pH 7.4, containing enzyme conjugate stabilizer, and 0.1% proclin as preservative. The protein concentration of peroxidase conjugated avidin is 2.5 mg/ml.

AVIDIN (from egg white) is conjugated with highly purified horseradish peroxidase (RZ > 3.0)

RECOMMENDED WORKING DILUTION FOR ELISA:

1:2,000 to 1:5,000. Working dilutions for specific applications should be determined by the investigator to obtain the best conditions.

STORAGE:

PI434423

Store at 2-8°C.

Explanation of symbols

Symbol	Description	Symbol	Description
REF	Catalogue Number	LOT	Batch code
RUO	Research Use Only	IVD	In vitro diagnostic medical device
\overline{X}	Use by	ł	Temperature limitation
***	Manufacturer	EC REP	European Community authorised representative
[-]	Without, does not contain	[+]	With, contains
Street Light	Protect from light	Æ	Consult accompanying documents
$\prod i$	Directs the user to consult instructions for use (IFU), accompanying the product.		

For research use only. CAUTION: Not for human or animal therapeutic or diagnostic use.

www.invitrogen.com

Invitrogen Corporation • 542 Flynn Rd • Camarillo • CA 93012 • Tel: 800.955.6288

E-mail: techsupport@invitrogen.com

(Rev 03/10) DCC-10-0838

Important Licensing Information - These products may be covered by one or more Limited Use Label Licenses (see the Invitrogen Catalog or our website, www.invitrogen.com). By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses. Unless otherwise indicated, these products are for research use only and are not intended for human or animal diagnostic, therapeutic or commercial use.