Technical Data Sheet

PE-CF594 Mouse Anti-NHP CD45

Product Information

Material Number:	562394
Alternate Name:	Pan Leukocyte, NHP-specific
Size:	50 tests
Vol. per Test:	5 µl
Clone:	D058-1283
Immunogen:	Rhesus peripheral whole blood
Isotype:	Mouse IgG1, κ
Reactivity:	QC Testing: Rhesus Macaque, Cynomolgus or Baboon
Storage Buffer:	Aqueous buffered solution containing BSA and $\leq 0.09\%$ sodium azide

Description

D058-1283 is a CD45 monoclonal antibody specific for non-human primate leukocytes, developed using Rhesus peripheral whole blood as immunogen. It does not cross react with human leukocytes. This antibody reacts with Baboon, Rhesus and Cynomolgus Macaque leukocytes in a similar pattern as seen with CD45 binding to the Leukocyte Common Antigen on human cells. Immunophenotypic analysis shows that D058-1283 binds to lymphocytes, monocytes and granulocytes of non-human primate blood samples. This antibody is able to block the binding of monoclonal antibody TÜ116; a reported anti-human CD45 that cross-reacts with non-human primate leukocytes. In Western blot analysis D058-1283 identifies a band of approximate molecular weight 180-200 kDa.

This antibody is conjugated to BD HorizonTM PE-CF594, which has been developed exclusively by BD Biosciences as a better alternative to PE-Texas Red[®]. PE-CF594 excites and emits at similar wavelengths to PE-Texas Red[®] yet exhibits improved brightness and spectral characteristics. Due to PE having maximal absorption peaks at 496 nm and 564 nm, PE-CF594 can be excited by the blue (488-nm), green (532-nm) and yellow-green (561-nm) lasers and can be detected with the same filter set as PE-Texas Red[®] (eg 610/20-nm filter).



Flow cytometric analysis of CD45 expression on Rhesus macaque peripheral blood leukocytes. Rhesus macaque whole blood was stained with BD Horizon[™] PE-CF594 Mouse Anti-NHP CD45 antibody (Cat. No. 562394; Right Panel) or with a BD Horizon[™] PE-CF594 Mouse IgG1, κ Isotype Control (Cat. No. 562292; Left Panel). The erythrocytes were lysed with BD Pharm Lyse[™] Lysing Buffer (Cat. No. 555899). Two parameter flow cytometric dot plots showing the correlated expression of CD45 (or Ig Isotype control staining) versus side-scattered light signals were derived from events with the forward light scatter signals of viable leukocytes. Flow cytometry was performed using a BD[™] LSR II Flow Cytometer System.

Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The antibody was conjugated with BD Horizon[™] PE-CF594 under optimum conditions, and unconjugated antibody and free PE-CF594 were removed.

BD	Bio	scier	nce

bdbiosciences.com				
United States Cana 877.232.8995 800.9	da Europe 979.9408 32.53.720.550	Japan 0120.8555.90	Asia Pacific 65.6861.0633	Latin America/Caribbean 55.11.5185.9995
For country contact in	nformation, visit bdbiosci	ences.com/conta	ct	
Conditions: The informati of any patents. BD Bioscie use of our products. Purch product or as a componer written authorization of I For Research Use Only. No	ion disclosed herein is not to b nnces will not be help responsi hase does not include or carry nt of another product. Any use Becton, Dickinson and Compar to ruse in diagnostic or ther.	e constructed as a re ble for patent infring any right to resell or e of this product oth ny is stictly prohibite apeutic procedures.	ecommendation to u gement or other vio r transfer this produc er than the permitte d. Not for resale.	ise the above product in violatio lations that may occur with the ct either as a stand-alone id use without the express

Application Notes

Application		
Flow cytometry	Routinely Tested	

Suggested Companion Products

Catalog Number	Name	Size	Clone
562292	PE-CF594 Mouse IgG1, κ Isotype Control	0.1 mg	X40
555899	Lysing Buffer	100 ml	(none)
554656	Stain Buffer (FBS)	500 ml	(none)

Product Notices

- This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1×10^{6} cells in a 100-µl experimental 1. sample (a test).
- 2. An isotype control should be used at the same concentration as the antibody of interest.
- 3. Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- 4. Please observe the following precautions: Absorption of visible light can significantly alter the energy transfer occurring in any tandem fluorochrome conjugate; therefore, we recommend that special precautions be taken (such as wrapping vials, tubes, or racks in aluminum foil) to prevent exposure of conjugated reagents, including cells stained with those reagents, to room illumination.
- 5. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- For fluorochrome spectra and suitable instrument settings, please refer to our Multicolor Flow Cytometry web page at 6. www.bdbiosciences.com/colors
- 7. CF[™] is a trademark of Biotium, Inc.
- 8. When excited by the yellow-green (561-nm) laser, the fluorescence may be brighter than when excited by the blue (488-nm) laser.
- 9 This product is provided under an Agreement between BIOTIUM and BD Biosciences. The manufacture, use, sale, offer for sale, or import of this product is subject to one or more patents or pending applications owned or licensed by Biotium, Inc. This product, and only in the amount purchased by buyer, may be used solely for buyer's own internal research, in a manner consistent with the accompanying product literature. No other right to use, sell or otherwise transfer (a) this product, or (b) its components is hereby granted expressly, by implication or by estoppel. This product is for research use only. Diagnostic uses require a separate license from Biotium, Inc. For information on purchasing a license to this product including for purposes other than research, contact Biotium, Inc., 3159 Corporate Place, Hayward, CA 94545, Tel: (510) 265-1027. Fax: (510) 265-1352. Email: btinfo@biotium.com.
- Because of the broad absorption spectrum of the tandem fluorochrome, extra care must be taken when using multi-laser cytometers, which 10. may directly excite both PE and CF™594.
- Texas Red is a registered trademark of Molecular Probes, Inc., Eugene, OR. 11
- Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols. 12.

References

Kishihara K, Penninger J, Wallace VA, et al. Normal B lymphocyte development but impaired T cell maturation in CD45-exon6 protein tyrosine phosphatase-deficient mice. Cell. 1993; 74(1):143-156. (Biology)

Knapp W, Dorken B, Rieber EP, et al, ed. Leucocyte Typing IV. New York: Oxford University Press; 1989:1-1208. (Biology)

Reimann KA, Waite BC, Lee-Parritz DE, et al. Use of human leukocyte-specific monoclonal antibodies for clinically immunophenotyping lymphocytes of rhesus monkeys. Cytometry. 1994; 17(1):102-108. (Biology)

BD Biosciences

babiosciences.	com	
United States	Canada	Europe
877.232.8995	800.979.9408	32.53.720.5

550 0120.8555.90 65.6861.0633 For country contact information, visit bdbiosciences.com/contact

Conditions: The information disclosed herein is not to be constructed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be help responsible for patent infringement or other violations that may occur with the use of our products. Buckness will not be help responsible for patent infringement or other violations that may occur with the product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Unless otherwise noted, BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2011 BD

Asia Pacific

Latin America/Caribbean

55.11.5185.9995

Japan

