Technical Data Sheet Jurkat + Pervanadate Lysate

Product Information		
Material Number:	611755	
Size:	500 μg	
Concentration:	1.0 mg/ml	
Storage Buffer:	SDS-PAGE buffer (62mM Tris pH 6.8, 2% SDS, 0.9% b-mercaptoethanol, 0.003% bromophenol blue, 5% glycerol)	

Description

Cell or tissue lysates are a useful control for the western blotting application. A source that can be rich in endogenously expressed proteins, this lysate is convenient to use and provided in a denatured format allowing for immediate loading onto an SDS-PAGE gel. Jurkat (T-cell leukemia) cells, were grown to 80-90% confluency, then treated with 1mM pervanadate, a general inhibitor of protein tyrosine phosphatases, for 15 minutes at 37°C.

Preparation and Storage

Store undiluted at -20° C. Do not dilute or boil.

Application Notes

Application

Western blot	Routinely Tested

Product Notices

- 1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- 2. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

BD Biosciences

bdbiosciences.com United States Canada Europe Japan 877.232.8995 888.259.0187 32.53.720.550 0120.8555.90 For country-specific contact information, visit bdbiosciences.com/how_to_order/

Asia Pacific 65.6861.0633

Latin America/Caribbean 55.11.5185.9995

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is stictly prohibited. For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2006 BD

