QUICK REFERENCE Transferrin Conjugates





Revision 1.0

Detailed protocol is available online at www.lifetechnologies.com/manuals.

Label	Cat. no.	Ex/Em maxima	Quantity
Biotin-XX	T23363	NA	5 mg
Fluorescein	T2871	494/518	5 mg
Alexa Fluor® 488	T13342	495/519	5 mg
Tetramethylrhodamine	T2872	555/580	5 mg
Alexa Fluor® 546	T23364	556/573	5 mg
Alexa Fluor® 555	T35352	555/565	5 mg
Alexa Fluor® 568	T23365	578/603	5 mg
pHrodo [™] Red	P35376	560/585	1 mg
Alexa Fluor® 594	T13343	590/617	5 mg
Texas Red®	T2875	595/615	5 mg
Alexa Fluor® 633	T23362	632/647	5 mg

For research use only. Not for use in diagnostic procedures.

WARNING! For every chemical, read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/support.

Quick Protocol

- 1. Plate cells and keep on ice for 10 minutes.
- Wash cells with cold Live Cell Imaging Solution (LCIS) (Cat. no. A14291DJ) containing 20 mM Glucose and 1% BSA.
- Add Transferrin conjugate at 25 µg/mL in LCIS containing 20 mM Glucose and 1% BSA. Incubate at 37°C for 5–20 minutes.
- Wash cells with LCIS (optional for the pHrodo[™] Red Transferrin conjugate).
- 5. Analyze cells using appropriate excitation/emission maxima given on page 1.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

DISCLAMER. LIFE TECHNOLOBIES CORPORATION AND/ORE ITS AFFLIATES) DISCLAM ALL WORAPATIES WITH RESPECT TO THIS DOCUMENT, DEPRESSED OR MAPLED, INCLUDING BUT NOT LIMITED TO THESE OF MERCHANTABILITY, FITNES FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT, TO THE EXTENT LIAUDED ST LAW, NOT DEVINT SHALL LIFE TENDICOGIES AND/OR ITS AFFLIARED ISE LABLE, WIETHER IN CONTRACT, TORY, WORAPATY, OR UNDER ANY STAUTE OR 0 ANY OTHER BASS FOR SPECIAL, INCOENTAL, INDIRECT, PUNTHE, MUTHEL OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARSING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOR

© 2013 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners.

For support visit lifetechnologies.com/support

lifetechnologies.com

7 October 2013