

Thermo Scientific Nalgene Quick Chill Unit

Cat. No. DS5114-0012

Application:

The Nalgene® Quick Chill™ Unit is designed to rapidly and repeatedly cool down biological samples such as ethanol precipitations for DNA/RNA, bacteria and viruses.

The self-contained unit eliminates the need to use traditional coolants such as dry ice/alcohol or acetone slurries which can be difficult to obtain and unsafe to handle.

There is no contact with the coolant so the label or marking stays intact. The unit can be used both at the lab bench and in the -70°C freezer.

On the bench at room temperature, the Quick Chill Unit will remain cold for 45 minutes for chilling multiple tubes.

See the next page for typical cooling times.

Adapters are included to convert any of the holes to accept 0.5-mL microcentrifuge tubes.

Thermo
SCIENTIFIC

Contact us for Sales and Service
[thermoscientific.com/contactus](https://www.thermoscientific.com/contactus)

*Contact information contained within
this document may be incorrect.

Thermo
SCIENTIFIC

© 2009 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

www.nalgenelabware.com

Asia: China Toll-free: 800-810-5118 or 400-650-5118;

India: +91 22 6716 2200, India Toll-free: 1 800 22 8374;

Japan: +81 3 3816 3355 Other Asian countries: 65 68729717

Europe: Austria: +43 1 801 40 0; Belgium: +32 2 482 30 30;

Denmark: +45 4631 2000; France: +33 2 2803 2180;

Germany: +49 6184 90 6940, Germany Toll-free: 08001-536 376;

Italy: +39 02 02 95059 or 434-254-375; Netherlands: +31 76 571 4440;

Nordic/Baltic countries: +358 9 329 100; Russia/CIS: +7 (812) 703 42 15;

Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 12;

UK/Ireland: +44 870 609 9203

North America: USA/Canada +1 585 586 8800;

USA Toll-free: 800 625 4327

South America: USA sales support: +1 585 899 7198

Countries not listed: +49 6184 90 6940 or

+33 2 2803 2180

8-0405-94 0609

Thermo
SCIENTIFIC