

Invitrogen™ Nucleic Acid Markers

Invitrogen offer the ultimate in convenience with both their Ready-to-Load Nucleic Acid Ladders.

The TrackIt™ and Ready-Load™ formats combine the most popular DNA ladders pre-mixed with loading buffer in a single tube. The TrackIt™ DNA ladders offer the added convenience of two tracking dyes, one that runs behind the sample and one that runs ahead of the sample allowing maximum resolution. Both formats offer:


- No thawing
- No diluting and mixing with loading dye
- No re-freezing

TrackIt™ format:

- Two tracking dyes which indicate when maximum resolution has been achieved
- Pre-mixed with loading buffer
- Ready to load
- Room temperature stability

Ready-Load™ format:

- Pre-mixed with loading buffer
- Ready to load
- Room temperature stability



* Stained with ethidium bromide