

Control

This product is prepared from purified plasma membranes isolated from an insect cell system (Sf9 cells infected with baculovirus), and can be used as negative control for other ABC transporter vesicle products.

* TMEP solution	
50 mM	Trls-HCl (pH 7.0)
50 mM	·······Mannitol
2 mM	·····EGTA
8 µg/ml	······Aprotinin
10 µg/ml ······	······Leupeptin
2 mM	DTT

Note

- Thaw vesicles rapidly at 37°C. Place on ice until use.
- It is recommended for reliable and reproducible data to minimize freeze/thaw cycle by dividing allquot to small portions.

THIS PRODUCT IS SUPPLIED FOR LABORATORY RESEARCH USE ONLY.

前田者

GenoMembrane, Inc. 75-1 Ono-cho, Tsurumi-ku, Yokohama, Kanagawa, 230-0046 Japan

TEL: +81-(0)45-508-2326, FAX: +81-(0)50-5548-3009,

URL: http://www.genomembrane.com, E-mail; info@genomembrane.com