

ELISA Wash Buffer - 10 x 1L Packets

Catalog Number: 00-0400

RUO: For Research Use Only

Product Information

Contents: ELISA Wash Buffer - 10 x 1L Packets

 Catalog Number: 00-0400



Temperature Limitation: Store at room temperature.



Batch Code: Refer to Vial



Use By: Refer to Vial

Description

The eBioscience ELISA Wash Buffer Powder is a dried powder formulation of phosphate buffered saline with 0.05% Tween-20. Each packet is reconstituted in 1 Liter of distilled, deionized water; the resulting 1x PBS/Tween-20 buffer (pH 7.4) is suitable for washing ELISA plate microwells. 10 packets of ELISA Wash Buffer Powder are supplied, yielding a total of 10 liters of 1x Wash Buffer.

Applications Reported

The eBioscience ELISA Wash Buffer Powder has been reported for use as the wash buffer for ELISA. 10 packets are supplied. Each packet is reconstituted in 1 liter of distilled, deionized water.

Applications Tested

The eBioscience ELISA Wash Buffer Powder has been tested by ELISA. The eBioscience ELISA Wash Buffer Powder is a dried powder formulation of phosphate buffered saline with 0.05% Tween-20. Each packet is reconstituted in 1 Liter of distilled, deionized water; the resulting 1x PBS/Tween-20 buffer (pH 7.4) is suitable for washing ELISA plate microwells. 10 packets of ELISA Wash Buffer Powder are supplied, yielding a total of 10 liters of 1x Wash Buffer. Detailed information and protocols about cytokine ELISA can be found in the BestProtocols® section.

Related Products

00-4201 1X TMB ELISA Substrate Solution

00-4202 ELISA Diluent Solution (5X)

00-4203 Super AquaBlue ELISA Substrate

18-4100 Avidin HRP

44-2404 Nunc MaxiSorp® flat-bottom 96 well plate

Not for further distribution without written consent.

Copyright © 2000-2010 eBioscience, Inc.

Tel: 888.999.1371 or 858.642.2058 • Fax: 858.642.2046 • www.eBioscience.com • info@eBioscience.com