

**3'-O-Me-m⁷G(5')ppp(5')G
RNA Cap
Structure Analog**



1-800-632-7799
info@neb.com
www.neb.com



S1411S 005120815081

S1411S

25 A₂₆₀ units Lot: 0051208

Store at -20°C Exp: 8/15

3'-O-Me-m⁷G(5')ppp(5')G Sodium Salt

Blocking of the 3'-hydroxyl of m⁷G with 3' O-Me assures that the capped transcripts are homogeneous. The 3'-hydroxyl of the non-methylated G is the only 3'-hydroxyl available for initiation (1).

Description: The 5' terminal m⁷G cap present on most eukaryotic mRNAs promotes translation in vitro at the initiation level (2,3,4). For most RNAs, elimination of the cap structure causes a loss of stability, especially against exonuclease degradation (5), and a decrease in the formation of the initiation complex of mRNAs for protein synthesis (5,6). Certain prokaryotic mRNAs containing a 5' terminal cap structure are translated as efficiently as or more efficiently than eukaryotic mRNAs in a eukaryotic cell-free protein synthesizing system (6). Also a cap requirement has been observed for splicing eukaryotic substrate RNAs (7).

A method using *E. coli* RNA Polymerase primed with m⁷G(5')ppp(5')G or m⁷G(5')ppp(5')A for an efficient in vitro synthesis of capped RNAs has been developed by Contreas (8). Larger amounts of capped RNAs are produced by transcription systems using SP6 RNA Polymerase primed with m⁷G(5')ppp(5')G (7).

Note: Addition of 131 µl water gives approximately a 10 millimolar solution.

Chromatographic Analysis:

HAISIL 300 C18 5 µm 50 x 10 mm
45 min grad linear 0.1 M TEAB 0–20% CH₃CN
RT = 10.9 minutes

TLC PEI Cellulose:

0.35 M LiCl 3.5 M urea
Mobility 0.77 vs xylene cyanol

Unit Definition:

MW = 817

ε₂₆₀ = ~19000

23.3 A₂₆₀ units/mg

25 A₂₆₀ units = ~1.07 mg = 1.31 micromoles and when dissolved in 131 µl water is approximately a 10 millimolar solution.

References:

1. Peng, Z.-H. et al. (2002) *ORG Lett.* 4(2)
2. Shatkin, A. J. (1978) *Cell.* 9, 645–653.
3. Fillipowicz, W. (1978) *FEBS Lett* 96, 1–11.
4. Banerjee, A.K. (1980) *Microbiol. Rev.* 44, 175–205.
5. Miura, K. (1981) *Adv. Biophys.* 14, 205–238.
6. Shatkin, A. J. et al. (1977) *Nucleic Acids. Res.* 4, 3065–3081.
7. Konarska, M. M. et al. (1984) *Cell* 38, 731–736.
8. Contreas, R. et al. (1982) *Nucleic Acids. Res.* 10, 6353–6363.
9. Paterson, B. M. and Rosenberg, M. (1979) *Nature* 279, 696–701.

CERTIFICATE OF ANALYSIS

**3'-O-Me-m⁷G(5')ppp(5')G
RNA Cap
Structure Analog**



1-800-632-7799
info@neb.com
www.neb.com



S1411S 005120815081

S1411S

25 A₂₆₀ units Lot: 0051208

Store at -20°C Exp: 8/15

3'-O-Me-m⁷G(5')ppp(5')G Sodium Salt

Blocking of the 3'-hydroxyl of m⁷G with 3' O-Me assures that the capped transcripts are homogeneous. The 3'-hydroxyl of the non-methylated G is the only 3'-hydroxyl available for initiation (1).

Description: The 5' terminal m⁷G cap present on most eukaryotic mRNAs promotes translation in vitro at the initiation level (2,3,4). For most RNAs, elimination of the cap structure causes a loss of stability, especially against exonuclease degradation (5), and a decrease in the formation of the initiation complex of mRNAs for protein synthesis (5,6). Certain prokaryotic mRNAs containing a 5' terminal cap structure are translated as efficiently as or more efficiently than eukaryotic mRNAs in a eukaryotic cell-free protein synthesizing system (6). Also a cap requirement has been observed for splicing eukaryotic substrate RNAs (7).

A method using *E. coli* RNA Polymerase primed with m⁷G(5')ppp(5')G or m⁷G(5')ppp(5')A for an efficient in vitro synthesis of capped RNAs has been developed by Contreas (8). Larger amounts of capped RNAs are produced by transcription systems using SP6 RNA Polymerase primed with m⁷G(5')ppp(5')G (7).

Note: Addition of 131 µl water gives approximately a 10 millimolar solution.

Chromatographic Analysis:

HAISIL 300 C18 5 µm 50 x 10 mm
45 min grad linear 0.1 M TEAB 0–20% CH₃CN
RT = 10.9 minutes

TLC PEI Cellulose:

0.35 M LiCl 3.5 M urea
Mobility 0.77 vs xylene cyanol

Unit Definition:

MW = 817

ε₂₆₀ = ~19000

23.3 A₂₆₀ units/mg

25 A₂₆₀ units = ~1.07 mg = 1.31 micromoles and when dissolved in 131 µl water is approximately a 10 millimolar solution.

References:

1. Peng, Z.-H. et al. (2002) *ORG Lett.* 4(2)
2. Shatkin, A. J. (1978) *Cell.* 9, 645–653.
3. Fillipowicz, W. (1978) *FEBS Lett* 96, 1–11.
4. Banerjee, A.K. (1980) *Microbiol. Rev.* 44, 175–205.
5. Miura, K. (1981) *Adv. Biophys.* 14, 205–238.
6. Shatkin, A. J. et al. (1977) *Nucleic Acids. Res.* 4, 3065–3081.
7. Konarska, M. M. et al. (1984) *Cell* 38, 731–736.
8. Contreas, R. et al. (1982) *Nucleic Acids. Res.* 10, 6353–6363.
9. Paterson, B. M. and Rosenberg, M. (1979) *Nature* 279, 696–701.

CERTIFICATE OF ANALYSIS