

# **Product Data Sheet**

## PE anti-human FOXP3 Flow Kit

Catalog # / Size: 320117 / 25 tests

Clone: 206D

**Isotype:** Mouse IgG1,  $\kappa$ 

Reactivity: Human, Cross-Reactivity: Baboon, Cynomolgus, Rhesus, Pigtailed

Macaque

Storage: This kit is guaranteed for three months. Upon receipt, store at 4°C and

protected from prolonged exposure to light. Do not freeze.

# **Applications:**

Applications: ICFC - Quality tested

**Recommended Usage:** Each lot of this antibody is quality control tested by immunofluorescent intracellular staining with flow cytometric analysis. For immunofluorescent

intracellular staining with flow cytometric analysis. For immunofluorescent staining, the suggested use of this reagent is ≤0.3 μg per million cells in 100 μl volume. It is recommended that the reagent be titrated for optimal

performance for each application.

**Application Notes: Materials Provided:** 

1. PE anti-human FOXP3, 25 tests

2. PE Mouse IgG1, k isotype control (ICFC), 25 tests

3. FOXP3 Fix/Perm buffer set, 100 tests (Cat. No. 421403)

Materials not included:

Cell Staining Buffer (Cat. No. 420201)

### **Surface Staining & FOXP3 Buffer Preparation:**

Centrifugation steps: perform at 250Xg for 5min Incubation steps: perform at room temperature

1. Perform cell surface staining if necessary (See protocol: Cell Surface Immunofluorescence Staining Protocol).

**2.** Prepare 1X buffer solutions of FOXP3 Fix/Perm buffer and FOXP3 Perm buffer in PBS.

**NOTE**: The FOXP3 Perm buffer (10X) may have crystalization or precipitation observed when it is stored at 2-8°C, however, it is normal and does not affect the buffer performance. If there is a heavy precipitation observed after diluting to 1X working solution, it may be clarified by filtering. **Caution:** The FOXP3 Fix/Perm buffer contains paraformaldehyde, which is

**Caution:** The FOXP3 Fix/Perm buffer contains paraformaldehyde, which is toxigenic and mutagenic. Please handle with caution and wear gloves, lab coat and necessary protection to avoid direct body contact.

### **FOXP3 Intracellular Staining Procedures:**

**3.** Add 1 ml of 1X Biolegend's FOXP3 Fix/Perm solution to each tube, resuspend the cells (gentle vortex) and incubate at room temperature in the dark for 20 minutes, then centrifuge and remove the supernatant. *The cell pellet will now be translucent and difficult to see; take care not to dislodge and accidentally aspirate cells at all later stages of staining protocol.* 

**4.** Wash: resuspend cells in cell staining buffer (Cat. No. 420201); centrifuge, then discard the supernatant.

**5.** Wash: resuspend in 1ml 1X BioLegend's FOXP3 Perm buffer; centrifuge, then discard the supernatant.

**6.** Resuspend cells in 1ml 1X BioLegend's FOXP3 Perm buffer, incubate in the dark for 15 minutes; centrifuge, then discard the supernatant. Resuspend the pellet in 100  $\mu$ l of 1X BioLegend's FOXP3 Perm buffer.

7. Add appropriate amount of flurochrome conjugated anti-FOXP3 antibody and incubate at room temperature in the dark for 30 minutes.


**8.** Wash twice with cell staining buffer (see step 4) then resuspend in 0.5 ml cell staining buffer. Analyze with flow cytometer using appropriate instrument settings.

NOTE: BioLegend's FOXP3 Fix/Perm buffer set (Cat. No. 421403) is



For research use only. Not for diagnostic use. Not for resale. BioLegend will not be held responsible for patent infringement or other violations that may occur with the use of our products.





turnan perpheral blood lymphocytes surface stained with CD4 FITC and then intracellular stained with 206D PE by using PE anti-human FOXP3 Flow Kit. Quadrant markers were set based on the staining of PE mlgG1, κ isotype control.

specifically developed and formulated for intracellular staining FOXP3 with minimum effect on surface fluorochrome staining and is highly recommended for optimal result of FOXP3 intracellular immunofluorescence staining.

#### **Application References:**

- 1. Roncador G, et al. 2005 Eur. J. Immunol. 35:1681. 2. Yang ZZ, et al. 2006. Blood 107:3639.
- 3. Liu W, et al. 2006. J. Exp. Med. 203:1701. PubMed Bollyky PL, et al. 2007. J. Immunol. 179:744.
   Bell MP, et al. 2007. J. Immunol. 179:1893.
   Tran DQ, et al. 2007. Blood doi:10.1182/blood-2007-06-094656. PubMed
- 8. Gao Q, et al. 2007. J Clin Oncol. 25: 2586. PubMed 9. Pillai V, et al. 2008. Blood 111:463. PubMed
- 10. Zheng Y, et al. 2008. J. Immunol. 181:1683. PubMed 11. Zonios DI, et al. 2008. Blood112:287. PubMed

- 12. Kavanagh B, et al. 2008. BloodPubMed
  13. Nevala WK, et al. 2009. Clin Cancer Res. 15:1931. PubMed 14. Grant J, et al. 2009. Cytometry B Clin Cytom. 76:69. PubMed 15. Nigam P, et al. 2010. J. Immunol. 184:1690. PubMed 16. Kmieciak M, et al. 2009. J. Transl. Med. 7:89. (ICFC) PubMed 17. Hartigan-O'Connor D. Jet al. 2007. J Exp Med. 204:2679. PubMed 2009. PubMed 2

- 18. Raghaven S, et al. 2009. Ann Rheum Dis. 68:1908. PubMed

Description: FOXP3 is a 50-55 kD transcription factor, also known as Forkhead box protein P3, Scurfin, JM2, or IPEX. It is proposed to be a master regulatory gene and more specific marker of T regulatory cells than most cell surface markers (such as CD4 and CD25). Transduced expression of FOXP3 in CD4 +/CD25 cells has been shown to induce GITR, CD103, and CTLA4 and impart a T regulatory cell phenotype. FOXP3 is mutated in X-linked autoimmunity-allergic dysregulation syndrome (XLAAD or IPEX) in humans and in "scurfy" mice. Overexpression of FOXP3 has been shown to lead to a hypoactive immune state suggesting that this transcriptional factor is a central regulator of T cell activity. In human, unlike in mouse, two isoforms of FOXP3 have been reported: one (FOXP3) corresponding to the canonical full-length sequence; the other (FOXP3  $\delta$ 2) lacking exon 2. The 206D antibody recognizes human FOXP3 epitope in the region of amino acids 105-235.

- Antigen References: 1. Hori S, et al. 2003. Science 299:1057.
  - 2. Gandhi R, et al. 2010. Nat. Immunol. 11:846.

Related Products: Product

Cell Staining Buffer Human TruŠtain FcX™ (Fc Receptor Blocking Solution) Clone Application FC, ICC, ICFC FC, ICC, ICFC



