

NativePAGE™ Running Buffer (20X)

(Novex®)

货号 BN2001

包装 1 L

NativePAGE™ Running Buffer (20X) is used to make the NativePAGE™ Cathode and Anode Running Buffers for use with an XCell SureLock® Mini Cell when running NativePAGE™ Gels.

详细说明

常用规格

Gel Type:	Native Gel
Buffer Type:	Running Buffers
Concentrated:	20 X
Product Size:	1 L
Gel Compatibility:	NativePAGE™ Gels
Shipping Condition:	Room Temperature
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

内容和储存

Store at room temperature