Novex® Tris-Glycine Native Sample Buffer (2X)

(Novex®)

货号 LC2673 包装 20 ml

描述

Novex® Tris-Glycine Native 2X Sample Buffer is ideal for native running conditions with NuPAGE® Tris-Acetate gels.

详细说明

常用规格

Gel Type:	Native Gel
Buffer Type:	Sample Loading Buffers
Concentrated:	2 X
Product Size:	20 ml
Gel Compatibility:	Novex®Tris-Glycine Gels

Shipping Condition: Wet Ice

Regulatory Statement:

For Research Use Only. Not for use in diagnostic procedures.

内容和储存

Store at 4°C to 25°C