

Novex® IEF Sample Buffer pH 3-10 (2X)

(Novex®)

货号 **LC5311**

包装 **25 ml**

描述

Novex® IEF Sample Buffer pH 3-10 (2X) is optimized pre-mixed IEF sample buffer for using with Novex® pH 3-10 IEF gels. IEF is a sensitive technique which is affected by many factors. The optimized pre-mixed IEF buffers reduce variability and ensure consistent results.

详细说明

常用规格

Gel Type:	IEF Gel
Buffer Type:	Sample Loading Buffers
Concentrated:	2 X
Product Size:	25 ml
Gel Compatibility:	Novex® IEF Gels
Shipping Condition:	Wet Ice
Regulatory Statement:	For Research Use Only. Not for use in diagnostic procedures.

内容和储存

Store in refrigerator (2–8°C).