

Human Erythropoietin Antibody

Polyclonal Rabbit IgG Catalog Number: AB-286-NA

DESCRIPTION

Species Reactivity	Human		
Specificity	Detects human Erythropoietin in direct ELISAs and Western blots. In direct ELISAs, approximately 35% cross-reactivity with recombinant rat Epo is observed, approximately 25% cross-reactivity with recombinant mouse Epo is observed, and approximately 10% cross-reactivity with recombinant canine Epo is observed.		
Source	Polyclonal Rabbit IgG		
Purification	Protein A or G purified		
Immunogen	Chinese hamster ovary cell line (CHO)-derived recombinant human Erythropoietin		
Endotoxin Level	<0.10 EU per 1 µg of the antibody by the LAL method.		
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details.		

APPLICATIONS

DATA

Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	Recombinant Human Erythropoietin (Ultrapure) (Catalog # 286-EP)
Immunohistochemistry	5-15 μg/mL	See Below
Neutralization	line. Kitamura, T. et	lity to neutralize Erythropoietin-induced proliferation in the TF-1 human erythroleukemic cell <i>al.</i> (1989) J. Cell Physiol. 140 :323. The Neutralization Dose (ND ₅₀) is typically 1-3 μg/mL in units/mL Recombinant Human Erythropoietin (Tissue Culture Grade).

Neutralization Immunohistochemistry Human Erythropoietin Antibody (µg/mL) Cell Proliferation Induced by Erythropoietin in Human Kidney Cancer Tissue. Erythropoietin and Erythropoietin was detected in immersion fixed paraffin-Neutralization by Human embedded sections of human kidney cancer tissue using 3000 3000 Erythropoietin Antibody. Rabbit Anti-Human Erythropoietin Polyclonal Antibody Recombinant Human (Catalog # AB-286-NA) at 10 µg/mL overnight at 4 °C 2500 2500 Erythropoietin (Tissue Culture Tissue was stained using the Anti-Rabbit HRP-DAB Cell & Grade) (Catalog # 287-TC) Tissue Staining Kit (brown; Catalog # CTS005) and Mean RFU Mean RFL stimulates proliferation in the counterstained with hematoxylin (blue). Lower panel shows a 2000 2000 TE-1 human erythroleukemic cell lack of labeling if primary antibodies are omitted and tissue is stained only with secondary antibody followed by incubation Protein Antiho line in a dose-dependent manner 1500 1500 with detection reagents. View our protocol for Chromogenic (orange line). Proliferation IHC Staining of Paraffin-embedded Tissue Sections. elicited by Recombinant Human 1000 1000 Erythropoietin (Tissue Culture 10-3 10-2 10-1 101 Grade) (0.3 units/mL) is neutralized (green line) by Recombinant Human Erythropoietin (units/mL) increasing concentrations of Rabbit Anti-Human Erythropoietin Polyclonal Antibody (Catalog # AB-286-NA). The ND₅₀ is typically 1-3 µg/mL

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 1 mg/mL in sterile PBS.		
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.		
Stability & Storage	ty & Storage Use a manual defrost freezer and avoid repeated freeze-thaw cycles.		
	 12 months from date of receipt, -20 to -70 °C as supplied. 		
	 1 month, 2 to 8 °C under sterile conditions after reconstitution. 		

• 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Erythropoietin (Epo) is a 34 kDa glycoprotein hormone in the type I cytokine family and is related to thrombopoietin. Its three N-glycosylation sites, four alpha helices, and N- to C-terminal disulfide bond are conserved across species. Glycosylation of Epo is required for biological activities *in vivo*. Mature human Epo shares 75%–84% amino acid sequence identity with bovine, canine, equine, feline, mouse, ovine, porcine, and rat EPO. Epo is primarily produced in the kidney by a population of fibroblast-like cortical interstitial cells adjacent to the proximal tubules. It is also produced in much lower, but functionally significant amounts by fetal hepatocytes and in adult liver and brain. Epo promotes erythrocyte formation by preventing the apoptosis of early erythroid precursors which express the Epo receptor (Epo R). Epo R has also been described in brain, retina, heart, skeletal muscle, kidney, endothelial cells, and a variety of tumor cells. Ligand induced dimerization of Epo R triggers JAK2-mediated signaling pathways followed by receptor/ligand endocytosis and degradation. Rapid regulation of circulating Epo allows tight control of erythrocyte production and hemoglobin concentrations. Anemia or other causes of low tissue oxygen tension induce Epo production by stabilizing the hypoxia-induceable transcription factors HIF-1α and HIF-2α. Epo additionally plays a tissue-protective role in ischemia by blocking apoptosis and inducing angiogenesis.

Rev. 5/2/2013 Page 1 of 1

