Phospho-PDGF Receptor α (Tyr762) Antibody

100 μl(10 western blots)

New 11/12

For Research Use Only. Not For Use In Diagnostic Procedures.

Applications W, IP Endogenous	Species Cross-Reactivity* H, (M, R)	Molecular Wt. 190 kDa	Source Rabbit**

Background: Platelet derived growth factor (PDGF) family proteins exist as several disulphide-bonded, dimeric isoforms (PDGF AA, PDGF AB, PDGF BB, PDGF CC, and PDGF DD) that bind in a specific pattern to two closely related receptor tyrosine kinases, PDGF receptor α (PDGFR α) and PDGF receptor β (PDGFR β). PDGFR α and PDGFR β share 75% to 85% sequence homology between their two intracellular kinase domains, while the kinase insert and carboxy-terminal tail regions display a lower level (27% to 28%) of homology (1). PDGFR α homodimers bind all PDGF isoforms except those containing PDGF D. PDGFRB homodimers bind PDGF BB and DD isoforms, as well as the PDGF AB heterodimer. The heteromeric PDGF receptor α/β binds PDGF B, C, and D homodimers, as well as the PDGF AB heterodimer (2). PDGFR α and PDGFR β can each form heterodimers with EGFR, which is also activated by PDGF (3). Various cells differ in the total number of receptors present and in the receptor subunit composition, which may account for responsive differences among cell types to PDGF binding (4). Ligand binding induces receptor dimerization and autophosphorylation. followed by binding and activation of cytoplasmic SH2 domain-containing signal transduction molecules, such as GRB2, Src, GAP, PI3 kinase, PLC_Y, and NCK. A number of different signaling pathways are initiated by activated PDGF receptors and lead

Provides PLCY, and NCK. A number of uniterint signaling pathways are initiated by activated PDGF receptors and lead to control of cell growth, actin reorganization, migration, and differentiation (5). Phosphorylation of PDGFR α at Tyr762 was identified at

Cell Signaling Technology using PTMScan[®], our LC-MS/ MS platform for phosphorylation site discovery (6). Tyr762 is located in the activation loop of the DDR1 kinase domain. Phosphorylation of PDGFR α at this site was also reported by several other labs to be a docking site for CrkII and CrkL upon induction by growth factor treatment (7,8).

Specificity/Sensitivity: Phospho-PDGF Receptor α (Tyr762) Antibody recognizes endogenous levels of PDGF receptor α protein only when phosphorylated at Tyr762. This antibody does not cross-react with endogenous levels of phosphorylated PDGF receptor β protein.



Western blot analysis of extracts from NCI-H1703 cells serum-starved overnight, untreated (-) or treated with Human Platelet-Derived Growth Factor AA (hPDGF-AA) #8913 (100 ng/ ml, 10 min;+), using Phospho-PDGF Receptor α (Tyr762) Antibody (upper) or PDGF Receptor α (D13C6) XP® Rabbit mAb #5241 (lower).

Source/Purification: Polyclonal antibodies are produced by immunizing animals with a synthetic phosphopeptide corresponding to residues surrounding Tyr762 of human PDGF receptor α protein. Antibodies are purified by protein A and peptide affinity chromatography.



Orders	877-616-CELL (2355)
	orders@cellsignal.com
Support	877-678-TECH (8324)
	info@cellsignal.com
Web	www.cellsignal.com

Entrez-Gene ID #5156 Swiss-Prot Acc. #P16234

Storage: Supplied in 10 mM sodium HEPES (pH 7.5), 150 mM NaCl, 100 μ g/ml BSA, 50% glycerol and less than 0.02% sodium azide. Store at -20°C. *Do not aliquot the antibody.*

*Species cross-reactivity is determined by western blot.

**Anti-rabbit secondary antibodies must be used to detect this antibody.

Recommended Antibody Dilutions: Western blotting

lestern blotting	 1:1000
nmunoprecipitation	1:100

For product specific protocols please see the web page for this product at www.cellsignal.com.

Please visit www.cellsignal.com for a complete listing of recommended complementary products.

Background References:

Ir

(1) Deuel, T.F. et al. (1988) *Biofactors* 1, 213-7.

(2) Bergsten, E. et al. (2001) Nat Cell Biol 3, 512-6.

(3) Betsholtz, C. et al. (2001) Bioessays 23, 494-507.

- (4) Coughlin, S.R. et al. (1988) Prog Clin Biol Res 266, 39-45.
- (5) Ostman, A. and Heldin, C.H. (2001) *Adv Cancer Res* 80, 1-38.
- (6) Rikova, K. et al. (2007) Cell 131, 1190-203.
- (7) Matsumoto, T. et al. (2000) *Biochem Biophys Res Commun* 270, 28-33.
- (8) Yokote, K. et al. (1998) Oncogene 16, 1229-39.

ц.

IMPORTANT: For western blots, incubate membrane with diluted antibody in 5% w/v BSA, 1X TBS, 0.1% Tween-20 at 4°C with gentle shaking, overnight.