Technical Data Sheet

PE Hamster Anti-Mouse CD61

Product Information

Material Number:	561910
Alternate Name:	Integrin β3 chain
Size:	25 μg
Concentration:	0.2 mg/ml
Clone:	2C9.G2
Immunogen:	Mouse T-cell Hybridoma 2B4 Vitronectin Receptor
Isotype:	Armenian Hamster IgG1, ĸ
Reactivity:	QC Testing: Mouse
	Tested in Development: Rat
Storage Buffer:	Aqueous buffered solution containing $\leq 0.09\%$ sodium azide.

Description

The 2C9.G2 antibody reacts with the integrin β 3 chain (CD61), which associates with the integrin α v chain (CD51) to form the vitronectin receptor and with the α IIb chain (CD41) to form the gpIIb/IIIa complex. Both receptors mediate adhesion to fibronectin, fibrinogen, vitronectin, thrombospondin, and von Willebrand factor. Leukocyte-endothelial adhesion is also mediated by the binding of $\alpha\nu\beta$ 3 integrin or vitronectin receptor to CD31 (PECAM-1). In addition, interaction of the $\alpha\nu\beta$ 3 integrin with its ligands regulates the L-type Ca2+ channel in vascular smooth muscle cells, possibly mediating vasodilatory responses to injury. Soluble and insoluble 2C9.G2 mAb mimics the effect of the natural ligands in smooth muscle cells from rat cremaster arterioles. Furthermore, osteopontin, also named Eta-1, is a cytokine that binds to $\alpha\nu\beta3$. CD61 is expressed on platelets, activated T lymphocytes, polymorphonuclear granulocytes, and blastocysts. Cross-reactivity of mAb 2C9.G2 to rat mast cells and platelets has been observed by flow cytometric analysis. mAb 2C9.G2 has been demonstrated to block binding of rat and mouse cells to fibronectin.

Preparation and Storage

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

The antibody was conjugated with R-PE under optimum conditions, and unconjugated antibody and free PE were removed.

Application Notes

An	nlic	atio	n

Flow cytometry		Routinely Tested		
Suggested Compa	nion Products			
Catalog Number	Name	Size	Clone	
553972	PE Hamster IgG1 κ Isotype Control	0.1 mg	A19-3	-
554656	Stain Buffer (FBS)	500 ml	(none)	

Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.

- 2. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- 3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 4. Although hamster immunoglobulin isotypes have not been well defined, BD Biosciences Pharmingen has grouped Armenian and Syrian hamster IgG monoclonal antibodies according to their reactivity with a panel of mouse anti-hamster IgG mAbs. A table of the hamster IgG groups, Reactivity of Mouse Anti-Hamster Ig mAbs, may be viewed at http://www.bdbiosciences.com/pharmingen/hamster chart 11x17.pdf.
- 5. For fluorochrome spectra and suitable instrument settings, please refer to our Fluorochrome Web Page at www.bdbiosciences.com/colors.
- 6. An isotype control should be used at the same concentration as the antibody of interest.

References

Ashkar S, Weber GF, Panoutsakopoulou V, et al. Eta-1 (osteopontin): an early component of type-1 (cell-mediated) immunity. Science. 2000; 287(5454):860-864. (Clone-specific: Blocking)

Frieser M, Hallmann R, Johansson S, Vestweber D, Goodman SL, Sorokin L. Mouse polymorphonuclear granulocyte binding to extracellular matrix molecules involves beta 1 integrins. *Eur J Immunol.* 1996; 26(12):3127-3136. (Biology)

Kieffer N, Phillips DR. Platelet membrane glycoproteins: functions in cellular interactions. Annu Rev Cell Biol. 1990; 6:329-357. (Biology)

Moulder K, Roberts K, Shevach EM, Coligan JE. The mouse vitronectin receptor is a T cell activation antigen. J Exp Med. 1991; 173(2):343-347. (Biology)

Piali L, Hammel P, Uherek C, et al. CD31/PECAM-1 is a ligand for alpha v beta 3 integrin involved in adhesion of leukocytes to endothelium. J Cell Biol. 1995; 130(2):451-460. (Clone-specific: Blocking)

BD Biosciences

bdbiosciences.c	om				
United States 877.232.8995	Canada 888.268.5430	Europe 32.53.720.550	Japan 0120.8555.90	Asia Pacific 65.6861.0633	Latin America/Caribbean 0800.771.7157
For country-spe	cific contact inf	formation, visit	bdbiosciences.co	m/how_to_orde	r/
Conditions: The in of any patents. BD use of our product product or as a co written authorizat For Research Use (BD, BD Logo and a	formation disclose Diosciences will not ts. Purchase does n mponent of anoth- tion of Becton Dick Dnly. Not for use in all other trademark	d herein is not to b ot be held responsi ot include or carry er product. Any us inson and Compan diagnostic or thera s are the property	e construed as a rec ble for patent infrin any right to resell ou e of this product oth y is strictly prohibite apeutic procedures. of Becton, Dickinsor	ommendation to use gement or other vio r transfer this produ her than the permitt ed. Not for resale. n and Company. ©20	e the above product in violation lations that may occur with the ct either as a stand-alone ed use without the express 111 BD

Schultz JF, Armant DR. Beta 1- and beta 3-class integrins mediate fibronectin binding activity at the surface of developing mouse peri-implantation blastocysts. Regulation by ligand-induced mobilization of stored receptor. *J Biol Chem.* 1995; 270(19):11522-11531. (Clone-specific: Blocking) Wu X, Mogford JE, Platts SH, Davis GE, Meininger GA, Davis MJ. Modulation of calcium current in arteriolar smooth muscle by alphav beta3 and alpha5 beta1 integrin ligands. *J Cell Biol.* 1998; 143(1):241-252. (Biology) Yasuda M, Hasunuma Y, Adachi H, et al. Expression and function of fibronectin binding integrins on rat mast cells. *Int Immunol.* 1995; 7(2):251-258.

(Clone-specific: Blocking)