


Histone Modifications



The Best Antibodies for Epigenetic Marks

Description	Cat. No.
acetyl-Histone H2A	07-376
acetyl-Histone H2A (K5)	07-290
acetyl-Histone H2A (K9)	07-289
phospho-Histone H2A (S129)	07-745
phospho-Histone H2A (S137)	09-018
phospho-Histone H2A/H4 (S1)	07-146
ubiquityl-Histone H2A	05-678

Description	Cat. No.
acetyl-Histone H2B	07-373
acetyl-Histone H2B (K5)	07-382
acetyl-Histone H2B (K12)	07-336
acetyl-Histone H2B (K15)	07-343
acetyl-Histone H2B (K16)	07-341
acetyl-Histone H2B (K120)	07-564
phospho-Histone H2B (S14)	07-191
ubiquityl-Histone H2B	05-1312

Description	Cat. No.
acetyl-Histone H3 (K9)	07-352
acetyl-Histone H3 (K14)	07-353
Anti-Histone H3.3	09-838
dimethyl-Histone H3 (K79)	04-835
trimethyl-Histone H3 (K4)	07-473
trimethyl-Histone H3 (K9)	07-442
trimethyl-Histone H3 (K27)	07-449
phospho-Histone H3 (S10)	06-570

Description	Cat. No.
acetyl-Histone H4	06-866
acetyl-Histone H4 (K5)	07-327
acetyl-Histone H4 (K8)	07-328
acetyl-Histone H4 (K12)	07-595
acetyl-Histone H4 (K16)	07-329
monomethyl-Histone H4 (K20)	04-735
dimethyl-Histone H4 (K20)	07-367
trimethyl-Histone H4 (K20)	07-463



Download an interactive version of this poster at:
www.millipore.com/histonemodapp



Merck Millipore is a division of **MERCK**

For more information on epigenetics antibodies, kits assays and proteins,
use your smart phone to scan the 2D barcode below,
or visit: www.millipore.com/epigenetics

