Technical Data Sheet Biotin Anti-Human CD137

| Product Information | |
|---------------------|--|
| Material Number: | 552533 |
| Alternate Name: | 4-1BB |
| Size: | 0.5 mg |
| Concentration: | 0.5 mg/ml |
| Clone: | 4B4-1 |
| Isotype: | Mouse IgG1 κ |
| Reactivity: | QC Testing: Human |
| Workshop: | VI C-7 |
| Storage Buffer: | Aqueous buffered solution containing $\leq 0.09\%$ sodium azide. |

Description

Reacts with 4-1BB, a 30 kDa glycoprotein expressed on activated T cells. 4-1BB can be detected on Con A or PHA-stimulated peripheral blood T cells (CD4+ and CD8+) and on CEM cells (human T-cell leukemia) following 2 day PMA and ionomycin stimulation, but not on resting T cells. 4-1BB is reported to be a participant in T-cell activation, being responsible for rescuing T cells from activation-induced apoptosis, upregulation of Th1-type, T-helper cells, downregulation of Th2-type cytokine production, and induction of cell adhesion to fibronectin.

Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. The antibody was conjugated with biotin under optimum conditions, and unreacted biotin was removed. Store undiluted at 4° C and protected from prolonged exposure to light. Do not freeze.

Application Notes

| Application | |
|-----------------|---------------------------|
| ELISA Detection | Tested During Development |

Recommended Assay Procedure:

ELISA Detection: Biotinylated 4B4-1 anti-human 4-1BB antibody is useful as a detection antibody for a sandwich ELISA which measures 4-1BB protein levels. Biotinylated 4B4-1 anti-human 4-1BB antibody can be paired with the purified h4-1BB-M127 anti-human 4-1BB as the capture antibody, with recombinant human 4-1BB:Fc fusion protein as the standard. This detection should be titrated between $0.5-1 \mu g/ml$ to determine its optimal concentration for ELISA detection. To obtain linear standard curves, doubling dilutions of recombinant human 4-1BB ranging from 2000 to 15 pg/ml are recommended for inclusion in each ELISA plate.

Note: This ELISA antibody pair shows no cross-reactivity with the following recombinant human cytokines: IL-1RI, IL-1RII, IL-1α, IL-1β, IL-2, IL-2sR-α, IL-3, IL-4, IL-4sR, IL-5, IL-6, IL-6R gp130, IL-7, IL-8, IL-10, IL-12 p40, IL-12 p70, IL-13, IL-15, TNF-α, TNF-β, sTNFRII, GM-CSF, TGFβ, IFN γ, CD40L, TRAIL, FasL.

Suggested Companion Products

| Catalog Number | Name | Size | Clone |
|----------------|--|------|-------|
| 552532 | Purified anti-human 4-1BB (Capture antibody) | | |

Product Notices

- 1. This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1 X 10e6 cells in a 100-µl experimental sample (a test).
- 2. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- 3. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.
- 4. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

BD Biosciences

| bdbiosciences.c | om | | | | |
|--|---|--|---|---|--|
| United States 877.232.8995 For country-spe | Canada 888.259.0187 cific contact infor | Europe 32.53.720.550 mation, visit bdbio | Japan 0120.8555.90 osciences.com/how | Asia Pacific 65.6861.0633 v_to_order/ | Latin America/Caribbear 55.11.5185.9995 |
| of any patents. BD use of our product product or as a con written authorizati For Research Use C | Biosciences will not b s. Purchase does not i nponent of another p on of Becton Dickinso Inly. Not for use in dia | e held responsible for nclude or carry any rig roduct. Any use of th on and Company is str gnostic or therapeuti | r patent infringement o ght to resell or transfer is product other than t | | may occur with the a stand-alone |

References

Kim YJ, Broxmeyer HE. Therapeutic potential of 4-1BB (CD137) as a regulator for effector CD8(+) T cells. J Hematother Stem Cell Res. 2001; 10(4):441-449.

Kim YJ, Broxmeyer HE. Therapeutic potential of 4-100 (co.iv) as a log-line to see a log-line to the set of the