# **Technical Data Sheet**

# PE Goat anti-Rat Ig Detector

#### **Product Information**

**Material Number:** 558551 Size: 25 tests 50 µl Vol. per Test: Isotype: Goat Ig

Storage Buffer: Aqueous buffered solution containing BSA and ≤0.09% sodium azide.

#### Description

The PE Goat anti-Rat Ig Detector is useful in verification of the conjugation of rat antibody of IgG and IgM isotypes to the BD CBA Functional Beads. Weaker reactivity may be seen with rat IgG2b isotype antibodies. It has minimal cross-reactivity with mouse, human, bovine or horse serum proteins.

## **Preparation and Storage**

The polyclonal antibody was purified from antiserum by negative adsorption and affinity chromatography.

The antibody was conjugated with R-PE under optimum conditions, and unconjugated antibody and free PE were removed.

Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze.

## **Application Notes**

#### Application

Flow cytometry	Routinely Tested	

## **Recommended Assay Procedure:**

This polyclonal antibody has been formulated at 50 µl/test and is useful in detection of rat antibody conjugated to the BD CBA Functional Beads. For specific instructions on the use of this antibody for the detection of antibody conjugated to BD CBA Functional Beads, refer to the

Conjugation Confirmation procedure in the BD CBA Functional Bead Conjugation Buffer Set Instruction Manual (Cat. No. 558556).

# **Suggested Companion Products**

Catalog Number	Name	Size	Clone
558556	Functional Bead Conjugation Buffer Set	15 reactions	(none)

#### **Product Notices**

- This conjugated product is sold under license to the following patents: US Patent No. 5,798,276.
- Source of all serum proteins is from USDA inspected abattoirs located in the United States.
- Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

# **BD Biosciences**

bdbiosciences.com

Asia Pacific Europe 877.232.8995 888.268.5430 32.53.720.550 0120.8555.90 65.6861.0633 0800.771.7157

For country-specific contact information, visit bdbiosciences.com/how\_to\_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation cof any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2011 BD

558551 Rev. 2 Page 1 of 1