

### **References for Product 12513**

1. Wehrman TS, von Degenfeld G, Krutzik PO, Nolan GP, Blau HM. (2006) Luminescent imaging of beta-galactosidase activity in living subjects using sequential reporter-enzyme luminescence. *Nat Methods*, 3, 295.
2. Yang X, Janatova J, Andrade JD. (2005) Homogeneous enzyme immunoassay modified for application to luminescence-based biosensors. *Anal Biochem*, 336, 102.
3. Masuda-Nishimura I, Fukuda S, Sano A, Kasai K, Tatsumi H. (2000) Development of a rapid positive/absent test for coliforms using sensitive bioluminescence assay. *Lett Appl Microbiol*, 30, 130.
4. Geiger R, Schneider E, Wallenfels K, Miska W. (1992) A new ultrasensitive bioluminogenic enzyme substrate for beta-galactosidase. *Biol Chem Hoppe Seyler*, 373, 1187.